Signum HiLite pre[®]

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(GB) Working Instructions

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signum[®]

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Heraeus Kulzer GmbH
Grüner Weg 11
63450 Hanau (Germany)

info.lab@heraeus.com
www.heraeus-kulzer.com

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12. Dokumentenhistorie

01.01.2008 Erstausgabe

1. Scope of application	2
1.1 General	2
1.2 Description and type of device	2
1.3 EC Declaration of Conformity	2
1.4 Relevant EC Directives	2
2. Notes about safe use	3
2.1 Explanation of symbols	3
2.2 Damage in transit	3
2.3 User's duties	3
2.4 Logbook	3
2.5 Safety instructions	3
3. Intended use	4
4. Scope of supply	4
5. Device description	5
6. Setting up, installation and putting into operation	6
7. Working with the Signum HiLite pre	6
7.1 Charging unit	6
7.2 Polymerisation	7
8. Cleaning	8
9. Device disposal	8
10. Technical data	9
11. Notes about malfunctions	10
12. Document history	10

1. Scope of application

1.1 General

Author f48618; HKG – FEGE

These Working Instructions apply to

Art. No.	Type of equipment	Issue
66034243	Signum HiLite pre with mobile charging unit	01.01.2008

1.2 Description and type of device

Description of the machine	Machine type	Serial no.
Pregelling unit for intermediate polymerisation of light-curing composite materials in dental laboratories.	Signum HiLite pre	2007-10-0005 ff.

1.3 EC Declaration of Conformity

We, Heraeus Kulzer GmbH, Grüner Weg 11, 63450 Hanau, Germany, hereby declare that in terms of its design, type and quality we have placed on the market the device described below, Signum HiLite pre, conforms to the relevant fundamental safety and health requirements of the EC Directives (see chapter 1.4).

If the machine undergoes a modification which has not been agreed with us, this declaration becomes invalid.

In the event of faults or damage please call our free Service Hotline

on 08004372522

or contact your local field sales representative.

In other countries please contact your appropriate Heraeus Kulzer GmbH subsidiary or your local dealer.

Signum HiLite pre is a registered trademark belonging to Heraeus Kulzer GmbH.

1.4 Relevant EC Directives

EC Low-voltage Directive 73/23EEC with its amendments

Electromagnetic Safety 89/336 with its amendments

EN 61010-1:2001, EN 61010-2-0010:2003, EN 61326:1997 + A1:1988 + A2:2001 + A3:2003

2. Notes about safe use

2.1 Explanation of symbols

	Word	Explanation
	Important !!	Safety-relevant chapter and sections of the Working Instructions.
	Note !!	Note in the Working Instructions about optimum use of the device.
	Device disposal	German Electrical and Electronic Equipment Act ('ElektroG')

2.2 Damage in transit

After taking receipt of the delivery please check the device for damage in transit and report any damage to the carrier within 24 hours. Under no circumstances should you use a damaged device.

2.3 User's duties

In addition to compliance with the manufacturer's legal regulations it is the user's duty to ensure compliance and implementation of statutory obligations at the place of work, e.g. the obligation to provide training and compliance with the health and safety at work act, as well as all other applicable regulations and laws.

2.4 Logbook

We recommend keeping a logbook. The logbook is used to record inspections and all important work performed (e.g. repairs, changes).

2.5 Safety instructions

With this laboratory equipment the safety instructions concerning the protection of individuals, handling the device and the items being processed is largely dependent on the conduct of the persons using the device.

Before putting the device into operation, read the Working Instructions carefully and observe the information in order to avoid errors and any damage or injury caused thereby, especially damage to health.

With regard to setting up and operation, the respective national laws, specifications and guidelines must be observed in addition to the information contained in these Working Instructions. Work on electronic parts of the device may only be performed by trained service companies when it is in a safe state (de-energised). Only approved genuine spare parts and accessories may be used. The use of other parts presents unknown risks and is prohibited.

The operational reliability and safety of the device is only guaranteed if the necessary inspections, maintenance and repair work is performed by Heraeus Kulzer GmbH service partners or by staff who have received training from the manufacturer.

Heraeus Kulzer GmbH accepts no liability whatsoever for any damage caused to the device in an incident as a result of improper repair work which has not been performed by Heraeus Kulzer GmbH service partners or by staff who have received training from the manufacturer, or if for an exchange of parts no genuine spare parts or accessories are used.

3. Intended use

The device Signum HiLite pre is designed to be used for setting light-curing composite materials in dental laboratories.

Important !!


Do not use the device for any of the activities described below.

- Do not use the device to heat or prepare food
- Do not treat or use easily combustible materials, liquids or gases
- Do not use on patients

4. Scope of supply

- 1 x handpiece (incl. lithium ion battery) with Working Instructions
- 1 x charging unit, power supply with international plug set, glare shield, clip

5. Device description

Charging unit with handle

- 1 – Light exit window
- 2 – Groove for attaching the glare shield
- 3 – Slide switch (intermittent / continuous)
- 4 – Belt clip
- 5 – Groove for attaching to the charging unit
- 6 – Articulated arm
- 7 – Charging unit
- 8 – Red status LED
- 9 – Orange status LED
- 10 – Green status LED

6. Setting up, installation and putting into operation

Place the device on a firm, non-combustible surface in such a way that it is stable. For dimensions and weight see chapter 10 (Tech. Data).

No combustible materials must be deposited or stored in the direct vicinity of the charging unit.

After unpacking, check the device and accessories to make sure everything is present and dispose of the packaging in an environmentally-friendly manner.

See chapter 4 (Scope of supply).

Important !!

Safety regulations for setting up the device. It is absolutely essential that you observe the following instructions:

- Before putting into operation check to make sure that the figures on the rating plate agree with the local power supply figures.
- Check the mains lead and plug for damage before use. If there is damage, the device must not be plugged into the mains.
- Do not short the electrical contacts on the handpiece with metal objects or liquids.
- Set up the device at a shockproof place which is protected against splashes of water and other liquids.
- Do not set up the device above or near sources of heat. When setting up, leave sufficient space round the device.
- Do not expose the device to direct sunlight or ultraviolet radiation!
- Do not set up the device in potentially explosive atmospheres where highly flammable air prevails or can arise (anaesthetic mixtures, oxygen, etc.).

In the event of faults or damage please call our Service Hotline, tel. 0 800-43 72 52 2 or contact your local field sales representative. In other countries please contact your appropriate Heraeus Kulzer GmbH subsidiary or your local dealer.

Plug the charging unit in at the mains in accordance with VDE regulations and the regulations issued by your local power supply company.

7. Working with the Signum HiLite pre

The Signum HiLite pre is supplied with power from a rechargeable lithium ion battery which is built into the handpiece.

7.1 Charging unit

Plug the power supply cable into the socket at the back of the charging unit (7) and plug in at the mains. The GREEN LED (10) lights up.

To charge the battery:

Place the handpiece in the receptacle of the charging unit (7). The orange LED (9) lights up. The charging phase is completed as soon as the orange LED (9) goes out. If an error occurs during the charging phase, the red LED (8) lights up (see chapter 11. Notes about malfunctions).

Important !!

Safety regulations for the battery

- Keep the battery out of children's reach.
- Only use genuine Heraeus Kulzer GmbH batteries.
- To charge the battery only use the Heraeus Kulzer GmbH charging unit (7). Do not charge the battery with universal battery chargers. Risk of explosion and fire!
- Dispose of the battery in a suitable manner in accordance with the laws applicable. Under no circumstances may it be disposed of along with domestic waste. The owner is liable for any damage caused by improper disposal of the battery.
- Do not use the battery for purposes other than those specified.

Risk of explosion !!

- Do not open, pierce or squash the battery; it contains toxic substances.
- Do not burn the battery or expose it to high temperatures: risk of explosion!
- Do not short the ends of the battery: risk of burning and fire!
- The battery of the Signum HiLite pre is not charged when it is delivered.

7.2 Polymerisation

Connect up the following accessories in order to use the Signum HiLite pre: Place the glare shield in the groove at the top (2) of the handpiece, see chapter 5 (Device description). For pregelling, the mobile handpiece can be removed from the charging unit (7) and guided freely over the material to be polymerised. The slide switch (3) to start light curing is located on the handpiece. Depending on which direction you move the slide switch in, you can choose whether light curing is intermittent or continuous.

With the aid of the articulated arms (6) on the charging unit (7) you can firmly position the handpiece above the material to be polymerised.

Note !!

For the curing times of intermediate/final polymerisation please refer to the separate working instructions for the particular material. For final polymerisation an additional device is required, e.g. Heraflash.

Important !!

Safety regulations for polymerisation

- Each time before using, always check to make sure that the device is working perfectly and that the accessories are operable. In the event of faults the device must not be used. Contact Customer Service if the device is affected by faults.
- Never point the light beam at anyone's eyes. Only point the light beam directly at the polymerisation material
- Before each curing cycle check to make sure that the glare shield is fitted to the top of the handpiece.
- When folding together, make sure you do not trap or squash your fingers.

Note !!

Avoid any direct contact between the Signum HiLite pre and the material to be polymerised. Polymerised composite residues on the light exit surface (1) impair light emission and have a negative effect on subsequent cycles.

8. Cleaning

Clean the unit with a moist cloth, making sure that no water gets into the housing of the device. Do not use aggressive detergents (abrasive powder, solvent).

Important !!

Safety regulations for cleaning

- The exterior of the charging unit (7) and the handpiece is not protected against the penetration of liquids! Avoid all contact between liquids and the surface and the electrical contacts!

9. Device disposal

German Electrical and Electronic Equipment Act ('ElektroG')

This law defines requirements for product responsibility in accordance with § 22 of the German Closed Substance Cycle and Waste Management Act for Electrical and Electronic Equipment. It chiefly aims at the avoidance of

waste electrical and electronic equipment and, in addition, at the re-use, recycling and other forms of utilisation of such wastes in order to reduce the quantity of waste to be disposed of and to reduce the emission of pollutants emerging from electrical and electronic equipment and wastes.

For disposal of used equipment of the type Signum HiLite pre please call our Service Hotline, tel. 0 800-43 72 52 2 or contact your local field sales representative.

In other countries please contact your appropriate Heraeus Kulzer GmbH subsidiary or your local dealer. The device must not be disposed of at local-authority waste disposal points.

10. Technical data

Power supply unit	
Rated input voltage	100–240 V AC, 0.6 A
Rated frequency	45–60 Hz
Rated output voltage	5 VDC
Charging unit	
Charging voltage	4.2 VDC, 1100 mA max.
Dimensions (W x H x D)	105 mm x 46 mm x 135 mm
Handpiece	
Power supply	Lithium ion battery, 3.7 V, 2200 mAh
Power LED	1.2 W
Wavelength	400 nm
Light curing	Intermittent or continuous
Battery charging time	3 hours, at nominal charging current
Dimensions and weight	L 140 mm x dia. 25 mm, 150 g

Ambient conditions:

- Temperature range 0°C to + 40°C
- Relative humidity between 30% and 75%
- Mains voltage $\pm$ 10% of nominal figure
- Use only indoors
- No direct ultraviolet radiation
- 2500 m max. above sea level
- Overvoltage category II
- Degree of pollution 2

	Fault	Possible cause	Remedy
Handpiece	No light curing	Battery not charged	Charge battery as described in chapter 7.1.
	Handpiece becomes warm	Continuous duty for 30 min	Switch off handpiece and allow to cool down
	Polymerisation is inadequate	Battery is down	Charge battery as described in chapter 7.1.
		Light exit window (1) soiled	Clean as described in chapter 8. Cleaning.

Base unit	Red LED (8) on charging unit (7) is lit	Overheating protection system activated	Allow the base unit to cool down
		Short circuit across the charging contacts	Battery faulty. Call the Service Hotline or contact your local field sales representative.
			Charging unit (7) faulty. Call the Service Hotline or contact your local field sales representative.
	The charging unit (7) cannot be switched on (green LED is not on)	The power supply cable is not connected up properly	Plug the cable into the charger (7) and the power socket
		Power supply cable is damaged	Change the power supply cable
		The charger (7) is not operating	Call the Service Hotline or contact your local field sales representative.

In the event of faults or damage please call our free Service Hotline, tel. 0 800-43 72 52 2 or contact your local field sales representative.
 In other countries please contact your appropriate Heraeus Kulzer GmbH subsidiary or your local dealer.

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